



- / We also offer blood collection needles for veterinary use. Please refer to our Website.
- / You can find more information about this product in our web-shop.



BLOOD COLLECTION NEEDLES

At Greiner Bio-One you will find a wide range of blood collection needles for human use but also for veterinary use.

Each product offers special properties and features that support you in your daily routine work. Blood collection needles

are intended for single use only. Our range includes Multiple Use Drawing Needles and Single Sample Needles as well as VACUETTE® VISIO PLUS Needles with optical puncture control. All cannulas are available in various dimensions and different lengths.

- / Sterile
- / Colour-coded needle dimensions
- / Available in different lengths
- / VISIO PLUS Needle with optical puncture control



VACUETTE® VISIO PLUS Needles

This cannula is equipped with a transparent plastic part in the middle, which serves as a viewing window. If the venipuncture is successful, this plastic part turns red due to the onset of blood flow.




Item no.	Colour code	Needle size	Needle length	Needle length	Qty. inner / outer
450040	● green	21 G	1.5 "	38 mm	2,000 / 6,000
450041	● black	22 G	1.5 "	38 mm	2,000 / 6,000
450042	● green	21 G	1 "	25 mm	2,000 / 6,000
450043	● black	22 G	1 "	25 mm	2,000 / 6,000



VACUETTE® Multiple Use Drawing Needles

Blood Collection Needles are predominantly used in routine blood collection and are exclusively single use. The colour of the needle protective cap simplifies visual recognition of the individual needle gauges.

Item no.	Colour code	Needle size	Needle length	Needle length	Qty. inner / outer
450071	● black	22 G	1 "	25 mm	2,000 / 6,000
450072	● green	21 G	1 "	25 mm	2,000 / 6,000
450073	● yellow	20 G	1 "	25 mm	2,000 / 6,000
450075	● black	22 G	1.5 "	38 mm	2,000 / 6,000
450076	● green	21 G	1.5 "	38 mm	2,000 / 6,000
450077	● yellow	20 G	1.5 "	38 mm	2,000 / 6,000
451071	● black	22 G	1 "	25 mm	100 / 1,000

Item no.	Colour code	Needle size	Needle length	Needle length	Qty. inner / outer
451072	 green	21 G	1 "	25 mm	100 / 1,000
451073	 yellow	20 G	1 "	25 mm	100 / 1,000
451075	 black	22 G	1.5 "	38 mm	100 / 6.000